



ELISA analysis of antibodies against vaccine-preventable pathogens



IgG antibodies against:

- Bordetella pertussis
- Diphtheria toxoid
- Measles virus
- Mumps virus
- Rubella virus
- TBE virus
- Tetanus toxoid
- Varicella zoster virus (VZV)

- Detailed standard curves for highest accuracy
- Reliable determination of low antibody titers
- Quantification in international units (rubella virus, measles virus, tetanus, diphtheria, B. pertussis, VZV) or Vienna units (TBE)
- Optimal monitoring of immune response using vaccination antigens (TBE, mumps virus)
- Fully automatable with EUROIMMUN Analyzer I and I-2P